ULTRASONIC CLEANER

VASCA ULTRASUONI

LINDA

3 e 6 L

Istruzioni d'uso / Instructions for use Mode d'emploi / Gebrauchsanleitung / Instrucciones de uso

INDICE

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TABLE OF CONTENTS

| | |
|--|---|
| CHAPTER 1 | 12 |
| Warranty | 12 |
| Services under warranty | 12 |
| Limitations on the warranty | 12 |
| Safety | 12 |
| Safety symbols | 12 |
| CHAPTER 2 | Errore. Il segnalibro non è definito. 13 |
| Electrical connection | 13 |
| Before operating the device | 13 |
| Chapter 3 | 15 |
| Command panel and description of symbols | 15 |
| Chapter 4 | 15 |
| Operating instructions | 15 |
| Chapter 5 | 16 |
| Accessories | 16 |
| Chapter 6 | 16 |
| Washing with ultrasonic | 16 |
| Useful tips and suggestions | 17 |
| Device maintenance | 17 |
| Chapter 7 | 18 |
| Assistance | 18 |
| Technical characteristics | 18 |
| Chapter 8 | 19 |
| Correct disposal of the product (electrical and electronic waste) | 19 |

Warranty

The information contained in this document is subject to possible variations without prior notice.

The device is guaranteed against defects regarding materials and manufacturing for a period of 24 months from the date of purchase according to the terms of the warranty certificate. During the warranty period, the manufacturer, in the event of device failure, will decide whether to repair or replace the defective product.

Services under warranty

For services or repairs under warranty, the device must be sent to FARO S.p.A via the dealer. The warranty does not cover the travel expenses of technicians, shipping costs and transport risks, which are borne by the buyer. In any case, the buyer will be responsible for all shipping costs (including taxes) of the products from another country.

Limitations on the warranty

The warranty covers the replacement or repair of defective components due to manufacturing defects, including the required labour. The above may not be applied to defects resulting from maintenance, improper or inappropriate use by the Customer, unauthorised modifications to the device, use of the product under environmental conditions other than those specified in this manual or by inadequate preparation of the installation site. No compensation will be due during the time required to replace or repair the device. Replacement will in any case be decided by the manufacturer and only if the device is verified to be totally unusable and repair not possible. Any direct or indirect damage of any nature to persons or property from the use or suspension of use of the device is excluded.

Safety

To preserve the original safety features of the product, the Customer must not replace any part or make any unauthorised changes.

WARNING

This symbol focuses attention on a procedure (s) whose failure or partial compliance may produce partial or total damage to the product or lead to injury to the operator. Before performing the procedures reported after this symbol, be sure to fully understand and comply with the specified conditions.

Safety symbols


Symbol of the instruction manual

Consult the manual to prevent any damage to the product or any physical damage to the operator.

High voltage symbol

Do not open or tamper with the device. Live electrical parts inside the device.

CHAPTER 2

Electrical connection

Before connecting the plug of the device to the power outlet, check that the voltage indicated on the nameplate corresponds to the voltage of your workplace. (See Fig. 2-1).

The earthing of the device is compulsory according to law. It is prohibited any intentional interruption of the earthing conductor inside or outside the device or the removal of the earth pin from the plug, as it makes the use of the device dangerous.

The manufacturer assumes no liability towards people or property resulting from failure to observe this rule. The electrical safety of this device is ensured only when it is properly connected to an effective earthing system, as required by the electrical safety regulations in force.

If you do not have an electrical system with regular earthing, do not connect the device to the power outlet and consult a qualified electrician as soon as possible.

WARNING

THIS DEVICE REQUIRES EARTHING.

Before operating the device

Before switching on the device fill the steel tank with water or cleaning fluid up to a maximum level of 3 cm from the top.

Fill the tank only if the power is disconnected. Accidental spillage of liquid may cause a short-circuit or electrocution and therefore be dangerous for the operator.

Make sure that the device is not damaged. Do not use damaged devices; if in doubt, ask your retailer or the manufacturer directly.

Install the ultrasonic cleaner on a flat and stable surface, able to support the weight of the device, utensils or objects that are inserted for the washing and the liquid. Handle it with care.

Install the ultrasonic washing machine away from heat sources. In addition, be careful not to install it near sources of moisture, wet work surfaces, or near sources of dust.

Make sure that the washing machine's feet are in the correct position to ensure air circulation.

Do not use the ultrasonic washing machine if: the power cord or plug is damaged, does not work properly, has been damaged, or has fallen. This may lead to electrical shock, fire or other accidents. Do not attempt to personally service the device. Contact only qualified technical personnel or the manufacturer.

Do not immerse the power cord or the plug in water. Keep the power cord away from hot surfaces.

Do not let the cord hang over the edges of tables or furniture.

Do not leave or do not use this device outdoors.

Do not tamper for any reason the electronic system of the ultrasonic washing machine "Electric shock hazard". For any repairs, always contact the manufacturer's technical service department.

Do not replace the power cord. If the cable is damaged due to


wear, or other reasons, turn off the device immediately and contact the service assistance department. **Do not lift and carry** the device when it is full of liquid. Handles are designed solely to lift and carry the device when it is empty and when disconnected from the power outlet.

WARNING

This ultrasonic washing machine works only with water or ultrasonic cleaning agents recommended by the manufacturer. Use only the amount of liquid needed to wash. **Do not use acidic or strong alkaline solutions** such as sodium hypochlorite because these substances, when in direct contact with the steel tray, cause a series of microscopic holes with irreparable and especially dangerous consequences for the operation of your device. As a matter of fact, all acidic or alkaline substances that can release corrosive compounds such as chlorine or other chemical compounds, or glutaraldehyde disinfectants, used together with the energetic ultrasonic cavitation activity, even cause accelerated corrosion of stainless steels. If it is necessary to use substances that could damage the steel tray, use glass or plastic cups. Remember that hydrogen fluoride breaks glasses and its use is extremely dangerous to health. Remember that any acid particles, micronised by intense ultrasonic cavitation activity, are dispersed in the working environment **causing corrosion to the device and also damaging health**. Do not use substances such as petrol, benzene, benzol or other harmful, explosive or flammable solvents. Use only solutions that are appropriate to the type of work to be performed. To prevent damage to the device: periodically replace the solution, do not operate the device without liquid, do not put parts or containers directly in contact with the bottom of the washing tray; use a rack or wire to keep the objects suspended. When the ultrasonic or heating is in operation, make sure that the liquid level does not drop by more than 1 cm below the operating level. Failure to follow such instructions may cause damage to the transducers or the heater, invalidating the warranty.

Important when using the ultrasonic washing machine

Before turning on the device fill the steel tray with water and/or detergent liquid to a maximum level of **3 cm** from the top; if your device version is equipped with the liquid drainage device, check that the top is properly closed. This will prevent accidental spillage of liquid substances.

Always pay attention to the liquid level so that it does not fall below 2/3 of the total height of the tray. If this occurs, the heating element, ultrasonic transducers and electronic circuitry may be seriously damaged.

Do not immerse your hands in the tray during operation.

Aggressive liquids that could damage the tray should be placed in a special glass which is then placed in the tray containing water that transmits ultrasonic waves to the glass and therefore to the solution it contains. Remember not to place the glass on the bottom of the tray: a special perforated support is available for such containers.

Use this device only for the professional use considered in this manual and for the purpose for which it was designed. This device is designed to wash and detach excess materials from instruments, prostheses, cutters, probes, tweezers, metal objects, metal parts in general, objects and glass tubes, crystals, stones and archaeological objects, rings, spectacles, parts difficult to access manually and so on.

In addition, LINDA washers offer ideal solutions for emulsion, accelerating chemical reactions, solution blends, degassing liquids, decontamination and sediment dissolution.

CHAPTER 3

Command panel and description of symbols

Fig. 1

LINDA device command panel

5. Switch/timer knob
6. MIN symbol (washing time expressed in minutes)
7. Symbol ∞ (infinite washing time)
8. Heating On/Off switch

Resistance symbol

Washing temperature symbol (degrees Celsius)

Fig. 2

LINDA device rear side view

3. Serial number plate and device identification
4. Power cord input

CHAPTER 4

Operating instructions

5. To turn ON the thermostated heating at 60°C, press the light switch indicated as Heater (Fig.1 - point 4)
6. Set the washing time by turning the Timer/Switch knob clockwise (Fig.1 - point 1)
7. At the end of the set time the machine will turn off automatically. To stop the operation of the machine at any time, simply turn the Timer/Switch knob to the "0" position.
8. To set a washing time of more than 15 minutes, turn the Timer/Switch knob counterclockwise to the " ∞ " position. This will allow you to set the washing time manually and turn the

appliance on and off according to user requirements. When the washing cycle is complete, to switch off the heating, press the Heater switch (Fig.1 - point 4).

CHAPTER 5

Accessories

Stainless steel cover

To be used to cover the ultrasonic tray. It is also useful if you start the heating cycle of the washing liquid. As a matter of fact, in order to reach the set washing temperature both efficiently and quickly, it is recommended to cover the tray with the cover. **Be very careful** with the condensation that is generated

as a result of the heating steam on the inside of the cover, because when opening, some liquid and condensation may pour into electrical parts such as the power outlet. It is therefore necessary to lift the cover vertically, so as to allow the condensation to drip only into the ultrasonic tray.

Rectangular stainless steel net rack

The stainless steel rack is a very useful accessory as it allows to optimise the washing of objects and to then drain them into the tray.

Glasses

Special glasses complete with adjustable black rubber support ring are available. They are used to save detergent liquid if you want to clean a small number of items, or if you want to wash with a detergent solution in a glass while simultaneously rinsing in the other glass. They also serve to wash objects with particularly aggressive solutions that can't not go directly into the steel tray.

Stainless steel support glass holder.

The stainless steel support glass holder is used to support the glasses during washing. Remember to adjust the glass support ring so that the bottom of the glass is always in contact with the liquid in the tray, but at the same time does not come in contact with the bottom of the steel tray.

CHAPTER 6

Washing with ultrasonic

Ultrasonic: what is it?

Ultrasonic is vibrations of a material medium similar to acoustic waves, but with frequencies that are too high to be heard by the human ear. The human hearing frequency limit varies from about 10 kHz to about 18 kHz. The audible frequency limit decreases with age. Specifically, the LINDA ultrasonic cleaners work at a frequency of 39 kHz, in other words a frequency far beyond that audible.

The LINDA ultrasonic tanks are made up of an ultrasonic generator and two or three piezoelectric transducers (depending on the model) that are placed on the external bottom of a stainless steel tray containing the washing liquid. The electronic generator produces a continuous signal at the frequency of 39 kHz and drives the piezoelectric transducers that transform the electrical signal into a mechanical vibration. This vibrational energy is transmitted to the tray liquid at a frequency of 39,000 oscillations per second. These pressure and depression oscillations create an enormous quantity of microbubbles inside the liquid that, by imploding in rapid succession, create immense impact energy between the cleaning fluid and the surface to be cleaned. This is called "cavitation" and enables an efficient and safe cleaning system with reduction in time.

Sweep System Technology

The LINDA models with Sweep System technology contain the best you can get from an ultrasonic washing machine. In fact, the brand new generator incorporates a special ultrasonic oscillator

with **Sweep System technology**. With this technology, the output frequency of the ultrasonic generator is modulated around a central frequency, so the transducers working at the frequency of 39 kHz are modulated at a frequency between 38 and 40 kHz. This frequency modulation offers the following advantages:

- **decreases** washing times;
- **prevents** the damage of delicate parts to be cleaned;
- **decreases** the effects of stationary waves increasing the distribution of ultrasonic energy in the washing liquid;
- **improves** cleaning results and facilitates the **cavitation** in liquids that are difficult to cavitate with traditional ultrasonic systems.

Useful tips and suggestions

The type of detergent to use. To wash any instrument or object, you should principally use tap water and dissolve in it a quantity of detergent or solution depending on the type of material residue to be removed from the object itself. Only use undiluted detergent if it is explicitly specified in the product's directions for use.

To wash many objects. Do not overfill the ultrasonic washing machine. Always place the objects in a logical manner at the bottom of the rectangular rack, place it in the tray and proceed with the appropriate washing. Too many items washed simultaneously reduces the effectiveness of cleaning with ultrasonic.

When objects are very dirty. Use long washing times with the heating of the washing liquid compatible with the object to be cleaned.

When to replace the cleaning solution. To achieve efficient cleaning, change the cleaning solution often, especially if you use it to clean surgical instruments.

Shape and size. There are no special recommendations for the shape of the objects to be cleaned, however avoid placing in the tank particularly heavy objects of very large dimensions. The liquid must always cover the entire object.

Washing time. The washing time always depends on the type of object, type and amount of material to be removed.

Device maintenance

Always disconnect the plug before carrying out any maintenance. To always ensure the safety of the device, all you need to do is regularly ensure that the enclosure and power cord of the device are not damaged.

Cleaning is the only maintenance normally required.

It must be carried out with the device disconnected electrically. Only use neutral detergents to clean the inside of the tray and a soft cloth to clean the outside such as the frame and the command panel.

Do not leave dirt deposits inside the tray.

If the liquid outlet hole is obstructed by dirt residues, clean with a non-pointed flexible stick, making sure that the rubber hose is not punctured.

Assistance

If the device does not work, perform the following checks:

- Check that the plug is properly plugged into the socket.
- Check the level of liquid in the tray.

WARNING

Services should only be carried out by technicians trained by the manufacturer. It is dangerous for any other person to carry out repairs. If assistance is required, contact your dealer as soon as possible and send the device directly to the manufacturer for repair, or contact directly:

FARO S.p.A.
Via Faro, 15
20876 - Ornago (MB)
Tel. +39 03968781
Fax +39 0396878515
Email: service@faro.it
<http://www.faro.it>

Technical characteristics

| Model | SONICA 2200 MHD S3 | SONICA 3200 MH S3 |
|-----------------------------|--|-------------------|
| | LINDA 3 L | LINDA 6 L |
| Power supply voltage | 230/240V - 50/60Hz | |
| Consumption | 130 W | 180 W |
| Consumption with heating | 305 W | 355 W |
| Weight in kg. (approx.) | 3.8 | 4.5 |
| External dimensions (in mm) | 270x170x260 | 325x270x260 |
| Tray dimensions (in mm) | 240x140x100 | 300x240x100 |
| Number of transducers | 2 | 3 |
| Environmental conditions | Temperature from 5 to 40°C; relative humidity 80% up to 31°C with linear decrease up to 50% at 40°C. | |
| Installation conditions | Class II according to EN 61010-1 | |

Information for use as "Medical Device"

Name: Ultrasonic Cleaning Device UMDNS-nomenclature (ECRI) 14-263

Field of application: washing and disinfection of surgical and dental instruments.

Classification (Medical Devices Directive 2007/47/EC, Appendix IX): Class I; active, non-invasive, non-implantable medical device.

Correct disposal of the product (electrical and electronic waste)

The symbol on the product or on the documentation indicates that the product must not be disposed of with other household waste at the end of the life cycle. To avoid any damage to the environment or health caused by inappropriate waste disposal, you are encouraged to separate this product from other types of waste and recycle it in a responsible manner to promote the sustainable reuse of material resources. The components of the device must be disposed of in a collection and recovery centre for plastic, metal and electronic components in order to be recycled.

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SOLTEC S.r.l.